

## **PRODUCT SPECIFICATION**

Name	Adrenocorticotropic Hormone Monoclonal Antibody (ACTH mAb)			
Code	HM141A			
RECEIVER INFORMATION				
Storage Temperature	2-8 °C (up to 3 Years)			
Storage Notes	Avoid freezing			
PRODUCT INFORMATION				
Description	Monoclonal antibody directed to ACTH			
Specificity	Specific for ACTH			
Information				
	Clone	Host /Sourc	ce Isotype	
	1141A	Mouse	lgG1	Affinity Chromatography
	1141B	Mouse	lgG1	Affinity Chromatography
Application				
	Pair	Capture	Detection	Recommended use
	Pair 1	1041A	1041B	CLIA/POCT
	Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.			
Presentation	Single homogenous batch, 10mM PBS, pH7.2, 0.09% Sodium Azide			
HEALTH AND SAFETY				
Warning	For in-vitro use and research or further processing only. This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1-1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.			