



PRODUCT SPECIFICATION

Name Anti-Mullerian Hormone Monoclonal Antibody

(AMH mAb)

Code HM188A

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years), -80 °C (up to 5

Years)

Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to AMH

Specificity Specific for AMH

Information

Clone	Host /Source	Isotype	Purification
1188A	Mouse	lgG2a	Protein G
1188B	Mouse	lgG1	Protein G

Application

Pair	Capture	Detection	Recommended use
Pair 1	1188A	1188B	CLIA

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but

use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

HEALTH AND SAFETY

Warning For in-vitro use and research or further processing only.

This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes

of water.