

PRODUCT SPECIFICATION

Name	D-Dimer Monoclonal Antibody (D-D mAb)					
Code	HM091A					
RECEIVER INFORMATION						
Storage Temperature	2-8 °C (up to 2 Years), below -20 °C (up to 3 Years)					
Storage Notes	Avoid repeat freeze thaw cycles					
PRODUCT INFORMATION						
Description	Monoclonal antibody directed to D-D					
Specificity	Specific for D-D					
Information						
	Clone	Host /Sou	urce Iso	otype	Purification	
	1091A	Mouse As	cites Ig	G2b	Protein G	
	1091B	Mouse As	cites lo	gG1	Protein G	
	1091C	Mouse As	cites lo	gG1	Protein G	
	1091E	Mouse As	cites lo	gG1	Protein G	
	1091F	Mouse Ascites		gG1	Protein G	
Application						
	Pair	Capture	Detection	Rec	ommended use	
	Pair 1	1091A	1091C		Biochemical	
	Pair 2	1091A	1091B		POCT	
	Pair 3	1091E	1		LETIA	
	Pair 3	1	1091F		LETIA	
	Each laboratory should determine an optimum working titer for use its particular application. Other applications have not been tested b use in such assays should not necessarily be excluded.					
Presentation	Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300					



Warning

For in-vitro use and research or further processing only. This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.