

PRODUCT SPECIFICATION

Name	Estradiol Monoclonal Antibody (E2 mAb)			
Code	HM177A			
RECEIVER INFORMATION				
Storage Temperature	2-8 °C (up to 2 Years), below -20 °C (up to 3 Years)			
Storage Notes	Avoid repeat freeze thaw cycles			
PRODUCT INFORMATION				
Description	Monoclonal antibody directed to E2			
Specificity	Specific for E2			
Information				
	Clone	Host /Source	lsotype	Purification
	1177A	Goat	lgG1	Protein G
Application				
	Pair	Capture D	etection F	Recommended use
	Pair 1	1177A		CLIA
	Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.			
Presentation	Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300			
HEALTH AND SAFETY				
Warning	For in-vitro use and research or further processing only. This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.			