



PRODUCT SPECIFICATION

Name Gastrin-17 Monoclonal Antibody (G-17 mAb)

Code HM013D

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 3 Years)

Storage Notes Avoid freezing

PRODUCT INFORMATION

Description Monoclonal antibody directed to G-17

Specificity Specific for G-17

Information

Clone	Host /Source	Isotype	Purification
1013E	Mouse	lgG2b	Protein A
1013F	Mouse	lgG1	Protein A

Application

Pair	Capture	Detection	Recommended use
Pair 1	1013E	1013F	POCT

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, PBS

HEALTH AND SAFETY

Warning For in-vitro use and research or further processing only.

Page: 1 of 1