



## PRODUCT SPECIFICATION

Name Helcobacter Pylori Monoclonal Antibody (HP mAb)

Code HM025D

RECEIVER INFORMATION

**Storage Temperature** 2-8 °C (up to 1 Week), below -20 °C (up to 3 Years)

Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

**Description** Monoclonal antibody directed to HP

**Specificity** Specific for HP

Information

Clone	Host /Source	Isotype	Purification
1025A	Mouse Ascites	lgG1	Protein A
1025B	Mouse Ascites	IgG1	Protein A

**Application** 

Pair	Capture	Detection	Recommended use
Pair 1	1025A	1025B	POCT

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Presentation** Single homogenous batch, 10mM PBS, pH7.4

**HEALTH AND SAFETY** 

**Warning** For in-vitro use and research or further processing only.

Page: 1 of 1