

## PRODUCT SPECIFICATION

Name Interleukin-6 Monoclonal Antibody (IL-6 mAb)

Code HM022D

RECEIVER INFORMATION

2-8 °C (up to 2 Years), below -20 °C (up to 3 Years), -80 °C (up to 5 **Storage Temperature** 

Years)

Avoid repeat freeze thaw cycles **Storage Notes** 

PRODUCT INFORMATION

Description Monoclonal antibody directed to IL-6

**Specificity** Specific for IL-6

Information

Clone	Host /Source	Isotype	Purification
1022F	Mouse	IgG1	Protein G
1022G	Mouse	IgG1	Protein G
1022H	Mouse	IgG1	Protein G
10221	Mouse	IgG2a	Protein G

**Application** 

Pair	Capture	Detection	Recommended use
Pair 1	1022F	1022G	CLIA
Pair 2	1022H	10221	POCT

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but

use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

**HEALTH AND SAFETY** 

Warning For in-vitro use and research or further processing only.

> This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes

of water.