

PRODUCT SPECIFICATION

Name	Lipoprotein-Assoeiated PhospholipaseA2 Monoclonal Antibody (Lp-PLA2 mAb)				
Code	HM003C				
RECEIVER INFORMATION					
Storage Temperature	2-8 °C (up to 2 Years), below -20 °C (up to 3 Years), -80 °C (up to 5 Years)				
Storage Notes	Avoid repeat freeze thaw cycles				
PRODUCT INFORMATION					
Description	Monoclonal antibody directed to Lp-PLA2				
Specificity	Specific for Lp-PLA2				
Information					
	Clone	Host /Source		lsotype	Purification
	1003A	Mouse		lgG2b	Protein G
	1003B	Mouse Mouse		lgG2b	Protein G
	1003C			lgG2b	Protein G
	1003D	Mous	e	lgG2b	Protein G
Application					
, pp. routien	Pair	Capture	Detection	Reco	mmended use
	Pair 1	1003A	1003C	Biochemical	
	Pair 2	1003B	1003D		POCT
	Each laboratory should determine an optimum working titer for use its particular application. Other applications have not been tested bu use in such assays should not necessarily be excluded.				
Presentation	Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300				
HEALTH AND SAFETY					
Warning	For in-vitro use and research or further processing only. This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.				
					Holmes Biotechnolog