

## **PRODUCT SPECIFICATION**

Name	Mycoplasma Pneumoniae Monoclonal Antibody (MP mAb)			
Code	HM183A			
RECEIVER INFORMATION				
Storage Temperature	2-8 °C (up to 1 Week), below -20 °C (up to 3 Years)			
Storage Notes	Avoid repeat freeze thaw cycles			
PRODUCT INFORMATION				
Description	Monoclonal antibody directed to MP			
Specificity	Specific for MP			
Information				
	Clone	Host /Sou	irce Isot	ype Purification
	1183A	Mouse Aso	cites IgC	61 Protein A
	1183B	Mouse Aso	cites IgC	G1 Protein A
Application				
	Pair	Capture	Detection	Recommended use
	Pair 1	1183A	1183B	CLIA/POCT
	Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.			
Presentation	Single homogenous batch, 10mM PBS, pH7.4			
HEALTH AND SAFETY				

Warning

For in-vitro use and research or further processing only.