

PRODUCT SPECIFICATION

Name	Procalcitonin Monoclonal Antibody (PCT mAb)			
Code	HM020B			
RECEIVER INFORMATION				
Storage Temperature	2-8 °C (up to 2 Years), below -20 °C (up to 3 Years)			
Storage Notes	Avoid repeat freeze thaw cycles			
PRODUCT INFORMATION				
Description	Monoclonal antibody directed to PCT			
Specificity	Specific for PCT			
Information				
	Clone	Host /Sour	rce Isoty	vpe Purification
	1020A	Mouse	lgG:	2a Protein G
	1020E	Mouse	lgG	1 Protein G
	10201	Mouse	lgG:	2b Protein G
Application				
	Pair	Capture	Detection	Recommended use
	Pair 1	1020A	10201	CLIA
	Pair 2	1020A	1020E	POCT
	Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.			
Presentation	Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300			
HEALTH AND SAFETY				
Warning	For in-vitro use and research or further processing only. This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.			