



## PRODUCT SPECIFICATION

Name Procollagen III N-Terminal Propeptide Monoclonal

**Antibody (PIIINP mAb)** 

Code HM027A

RECEIVER INFORMATION

**Storage Temperature** 2-8 °C (up to 1 Month), below -20 °C (up to 5 Years)

Storage Notes Avoid repeat freeze thaw cycles

**PRODUCT INFORMATION** 

**Description** Monoclonal antibody directed to PIIINP

**Specificity** Specific for PIIINP

Information

Clone	Host /Source	Isotype	Purification
1027A	Mouse	lgG1	Protein A
1027B	Mouse	lgG1	Protein A

**Application** 

Pair	Capture	Detection	Recommended use
Pair 1	1027A	1027B	CLIA/Biochemical

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Presentation** Single homogenous batch, 10mM PBS, pH7.4

**HEALTH AND SAFETY** 

**Warning** For in-vitro use and research or further processing only.

Page: 1 of 1