



## PRODUCT SPECIFICATION

Name Testosterone Monoclonal Antibody (T mAb)

Code HM115A

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years), -80 °C (up to 5

Years)

Storage Notes Avoid repeat freeze thaw cycles

**PRODUCT INFORMATION** 

**Description** Monoclonal antibody directed to T

**Specificity** Specific for T

Information

CloneHost /SourceIsotypePurification1115AGoatIgG1Protein G

**Application** 

Pair	Capture	Detection	Recommended use
Pair 1	1115A		CLIA

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but

use in such assays should not necessarily be excluded.

**Presentation** Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

**HEALTH AND SAFETY** 

**Warning** For in-vitro use and research or further processing only.

This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes

of water.

Page: 1 of 1