

## **PRODUCT SPECIFICATION**

Name	Thyroxine Monoclonal Antibody (T4 mAb)				
Code	HM151A				
RECEIVER INFORMATION					
Storage Temperature	2-8 °C (up to 2 Years), below -20 °C (up to 3 Years), -80 °C (up to 5 Years)				
Storage Notes	Avoid repeat freeze thaw cycles				
PRODUCT INFORMATION					
Description	Monoclonal antibody directed to T4				
Specificity	Specific for T4				
Information				-	
	Clone	Host /So	ource	Isotype	Purification
	1151A	Goa	t	lgG1	Protein G
Application					
	Pair	Capture	Detection	Recommended use	
	Pair 1	11	51A	CLIA	
Presentation	Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300				
Presentation	Single nomo	jenous batch,	20111111 PD3	, pπ <i>r</i> .2, υ.	
HEALTH AND SAFETY					
Warning	For in-vitro use and research or further processing only. This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.				