

PRODUCT SPECIFICATION

Name Cardiac Troponin I Monoclonal Antibody (cTnl mAb)

Code HM033E

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years), -80 °C (up to 5

Years)

Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to cTnI

Specificity Specific for cTnI

Information

Clone	Host /Source	Isotype	Purification
1033D	Mouse	IgG1	Protein G
1033E	Mouse	IgG1	Protein G
1033F	Mouse	lgG2b	Protein G

Application

Pair	Capture	Detection	Recommended use
Pair 1	1033D+1033E	1033F	CLIA

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but

use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

HEALTH AND SAFETY

Warning For in-vitro use and research or further processing only.

This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes

of water.