

PRODUCT SPECIFICATION

Name Soluble Suppression of Tumorigenicity 2

Monoclonal Antibody(sST2 mAb)

Code HM148A

RECEIVER INFORMATION

Storage Temperature 2-8 ° C (up to 2 Years), below -20 ° C (up to 3 Years), -80 ° C (up to

5 Years)

Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to sST2

Specificity Specific for sST2

Information

Clone	Host /Source	Isotype	Purification
1148A	Mouse	lgG1	Protein G
1148B	Mouse	lgG1	Protein G
1148C	Mouse	lgG1	Protein G
1148D	Mouse	lgG1	Protein G
1148E	Mouse	lgG1	Protein G
1148F	Mouse	lgG1	Protein G


Application

Pair	Capture	Detection	Recommended use
Pair 1	1148C	1148A	CLIA
Pair 2	1148B	1148F	POCT
Pair 3	1148[D/1148E	Others

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

Page: 1 of 2





HEALTH AND SAFETY

Warning

For in-vitro use and research or further processing only. This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.

Page: 2 of 2